

alpha Tubulin 8 antibody [N1C1-2]

Cat. No. GTX110432

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 50 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human alpha Tubulin 8. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX110432 IHC-P Image

Immunohistochemical analysis of paraffin-embedded SW480 xenograft, using alpha Tubulin 8(GTX110432) antibody at 1:500 dilution.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX110432 WB Image

Sample (30 ug of whole cell lysate)


A: NIH-3T3

10% SDS PAGE

GTX110432 diluted at 1:5000


GTX110432 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed A431, using alpha Tubulin 8(GTX110432) antibody at 1:200 dilution.


GTX110432 WB Image

Sample (30 ug of whole cell lysate)

A: A431 (GTX27909)

B: HeLa

10% SDS PAGE

GTX110432 diluted at 1:1000



For full product information, images and publications, please visit our [website](http://www.genetex.com).