

## XRCC2 antibody

Cat. No. GTX11052

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB
Reactivity	Human

References ( 1 )

Package

100 µg

## Applications

## Application Note

WB: Use at an assay dependent dilution. Predicted molecular weight: 33 kDa. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 32 kDa. ( [Note](#) )

Product Note Recognises as little as 10 ng of purified protein. Does not cross-react with Rad51B, Rad51C, Rad51D or XRCC3 by western.

## Properties

Form	Liquid
Buffer	PBS
Preservative	No preservative
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full length protein (Human).
Purification	Protein A purified
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).