

Human Caspase 3 protein, GST tag

Cat. No. GTX110543-pro
Species Human

Package

100 µg, 25 µg

PROPERTIES

Form Liquid

Buffer PBS, 0.1% SDS

Preservative No Preservative

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Region/Sequence N-terminal GST tagged partial fragment of NP_004337


Expression System E. coli

Purity >95%

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES



GTX110543-pro Image

10µg of GTX110543-pro Caspase 3 recombinant protein analyzed using SDS-PAGE and stained with coomassie blue and captured by black and white camera.


For full product information, images and publications, please visit our [website](#).