

Human HSBP1 protein, His tag

Cat. No. GTX110661-pro

Species	Human
----------------	-------

Package 100 µg, 25 µg

Properties

Form	Liquid
Buffer	PBS, 0.1% SDS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Region/Sequence	N-terminal His tagged full length protein of NP_001528
Expression System	E. coli
Purity	>90%
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).