

ARA9 antibody [N1C3]

Cat. No. GTX110665

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:10000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 38 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human ARA9. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx110665 WB Image

ARA9 antibody [N1C3] detects ARA9 protein by Western blot analysis.


A. 30 μ g HeLa whole cell lysate/extract

B. 30 μ g Molt-4 whole cell lysate/extract

C. 30 μ g Raji whole cell lysate/extract

10 % SDS-PAGE

ARA9 antibody [N1C3] (GTx110665) dilution: 1:5000




GTx110665 WB Image

Sample (50 μ g of whole cell lysate)

A: Mouse brain

10% SDS PAGE

GTx110665 diluted at 1:1000



GTx110665 WB Image

Various whole cell extracts (30 μ g) were separated by 10% SDS-PAGE, and the membrane was blotted with ARA9 antibody [N1C3] (GTx110665) diluted at 1:500.



For full product information, images and publications, please visit our [website](https://www.genetex.com).