

Human Aspartoacylase protein, His tag


Cat. No. GTX110699-pro

Species Human

References (1) Package 100 µg, 25 µg

Properties	
Form	Liquid
Buffer	PBS, 0.1% SDS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Region/Sequence	N-terminal His tagged full length protein of NP_001121557
Expression System	E. coli
Purity	>90%
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES



GTX110699-pro Image

 $10\mu g$ of GTX110699-pro Aspartoacylase recombinant protein analyzed using SDS-PAGE and stained with coomassie blue and captured by black and white camera.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 10 / 25 Page 1 of 1