

FBXL3 antibody

Cat. No. GTX110755

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human, Mouse

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000

Not tested in other applications.

Calculated MW 49 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.41 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human FBXL3. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx110755 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed A431, using FBXL3 (GTx110755) antibody at 1:200 dilution.




GTx110755 WB Image

FBXL3 antibody detects FBXL3 protein by Western blot analysis.
A. 30 µg C2C12 whole cell lysate/extract
10 % SDS-PAGE
FBXL3 antibody (GTx110755) dilution: 1:500



GTx110755 WB Image

Sample (30 ug of whole cell lysate)
A: Molt-4 (GTx27912)
10% SDS PAGE
GTx110755 diluted at 1:1000



GTx110755 WB Image

The data was published in the journal Nat Commun in 2017. [PMID: 28090088](https://pubmed.ncbi.nlm.nih.gov/28090088/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).