

GTPBP4 antibody

Cat. No. GTX110826

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human, Mouse

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
ICC/IF	1:100-1:1000

Not tested in other applications.

Calculated MW 74 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human GTPBP4. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2

DATA IMAGES

GTX110826 WB Image

GTPBP4 antibody detects GTPBP4 protein by western blot analysis.

A. 30 µg 293T whole cell lysate/extract


B. 30 µg A431 whole cell lysate/extract

C. 30 µg HeLa whole cell lysate/extract

D. 30 µg HepG2 whole cell lysate/extract

7.5 % SDS-PAGE

GTPBP4 antibody (GTX110826) dilution: 1:5000



GTX110826 ICC/IF Image

GTPBP4 antibody detects GTPBP4 protein at nucleolus by immunofluorescent analysis.

Sample: A431 cells were fixed in ice-cold MeOH for 5 min.

Green: GTPBP4 protein stained by GTPBP4 antibody (GTX110826) diluted at 1:500.

Blue: Hoechst 33343 staining.


GTX110826 WB Image

Sample (10 ug of whole cell lysate)

A: NIH-3T3

B: JC

C: BCL-1

7.5% SDS PAGE

GTX110826 diluted at 1:5000



For full product information, images and publications, please visit our [website](https://www.genetex.com).