

RPRM antibody [N1C3]

Cat. No. GTX110976

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-Fr
Reactivity	Human, Mouse

References (3)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-Fr	Assay dependent

Not tested in other applications.

Calculated MW 12 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.71 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human RPRM Recombinant protein.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX110976 WB Image

Non-transfected (–) and transfected (+) H1299 whole cell extract (30 µg) were separated by 12% SDS-PAGE, and the membrane was blotted with RPRM antibody [N1C3] (GTX110976) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.



GTX110976 WB Image

Non-transfected (–) and transfected (+) unboiled F1299 whole cell extracts (30 µg) were separated by 12% SDS-PAGE, and the membrane was blotted with RPRM antibody [N1C3] (GTX110976) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.



GTX110976 WB Image

The data was published in the journal PLoS One in 2015. [PMID: 26176219](https://pubmed.ncbi.nlm.nih.gov/26176219/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).