

Human DAPL1 protein, His tag

Cat. No. GTX111012-pro

Species	Human
---------	-------


Package

100 µg, 25 µg

Properties

Form	Liquid
Buffer	PBS, 0.1% SDS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Region/Sequence	N-terminal His tagged full length protein of NP_001017920
Expression System	E. coli
Purity	>95%
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES



GTX111012-pro Image

10µg of GTX111012-pro DAPL1 recombinant protein analyzed using SDS-PAGE and stained with coomassie blue and captured by black and white camera.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 1