

## Human NDUFA5 protein, His tag

Cat. No. GTX111016-pro

|         |       |
|---------|-------|
| Species | Human |
|---------|-------|

Package

100 µg, 25 µg

## Properties

|      |        |
|------|--------|
| Form | Liquid |
|------|--------|

|        |               |
|--------|---------------|
| Buffer | PBS, 0.1% SDS |
|--------|---------------|

|              |                  |
|--------------|------------------|
| Preservative | No preservatives |
|--------------|------------------|

|         |                                                                                                                                                   |
|---------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
|---------|---------------------------------------------------------------------------------------------------------------------------------------------------|

|               |                                                                                  |
|---------------|----------------------------------------------------------------------------------|
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.) |
|---------------|----------------------------------------------------------------------------------|

|                 |                                                        |
|-----------------|--------------------------------------------------------|
| Region/Sequence | N-terminal His tagged full length protein of NP_004991 |
|-----------------|--------------------------------------------------------|


|                   |         |
|-------------------|---------|
| Expression System | E. coli |
|-------------------|---------|

|        |      |
|--------|------|
| Purity | >95% |
|--------|------|

|             |              |
|-------------|--------------|
| Conjugation | Unconjugated |
|-------------|--------------|

|      |                                                                                                                                                   |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Note | For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------|

## DATA IMAGES



## GTX111016-pro Image

10µg of GTX111016-pro NDUFA5 recombinant protein analyzed using SDS-PAGE and stained with coomassie blue and captured by black and white camera.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 1