

TID1 antibody [N3C3]

Cat. No. GTX111077

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

References (3)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 52 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human TID1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX111077 IHC-P Image

Immunohistochemical analysis of paraffin-embedded MCF7 xenograft, using TID1(GTX111077) antibody at 1:500 dilution.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




GTX111077 IHC-P Image

TID1 antibody [N3C3] detects TID1 protein at cytoplasm in rat muscle by immunohistochemical analysis.

Sample: Paraffin-embedded rat muscle.


TID1 antibody [N3C3] (GTX111077) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX111077 ICC/IF Image

Immunofluorescence analysis of methanol-fixed MCF-7, using TID1(GTX111077) antibody at 1:500 dilution.



GTX111077 WB Image

Sample (30 ug of whole cell lysate)

A: 293T

B: A431 (GTX27909)

C: H1299

D: HeLa

E: Hep G2 (GTX27900)

F: Molt-4 (GTX27912)

G: Raji

10% SDS PAGE

GTX111077 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).