

EEF1D antibody

Cat. No. GTX111151

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 31 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human EEF1D. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX111151 IHC-P Image

Immunohistochemical analysis of paraffin-embedded human adenocarcinoma, using EEF1D(GTX111151) antibody at 1:250 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTX111151 WB Image

Sample (30 ug of whole cell lysate)

A: A549

10% SDS PAGE

GTX111151 diluted at 1:1000


GTX111151 WB Image

EEF1D antibody detects EEF1D protein by Western blot analysis.


A. 30 µg THP-1 whole cell lysate/extract

B. 30 µg HL-60 whole cell lysate/extract

C. 30 µg NCI-H929 whole cell lysate/extract

10 % SDS-PAGE

EEF1D antibody (GTX111151) dilution: 1:1000


GTX111151 ICC/IF Image

Immunofluorescence analysis of methanol-fixed HeLa, using EEF1D(GTX111151) antibody at 1:500 dilution.



For full product information, images and publications, please visit our [website](http://www.genetex.com).