

CIAPIN1 antibody

Cat. No. GTX111212

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 34 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human CIAPIN1 Recombinant protein.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTXT111212 ICC/IF Image


CIAPIN1 antibody detects CIAPIN1 protein at nucleus and cytoplasm by immunofluorescent analysis.

Sample: A375 cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: CIAPIN1 protein stained by CIAPIN1 antibody (GTXT111212) diluted at 1:500.

Blue: Hoechst 33342 staining.

Scale bar = 10 μ m.



GTXT111212 IHC-P Image

CIAPIN1 antibody detects CIAPIN1 protein at cytosol and nucleus on SkHep1 xenograft by immunohistochemical analysis.

Sample: Paraffin-embedded SkHep1 xenograft.

CIAPIN1 antibody (GTXT111212) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTXT111212 WB Image

Sample (30 ug of whole cell lysate)

A: Jurkat

B: Raji

C: K562


D: THP-1

E: HL-60

F: NCI-H929

10% SDS PAGE

GTXT111212 diluted at 1:1000


GTXT111212 WB Image

Sample (30 ug of whole cell lysate)

A: NIH-3T3

B: JC

C: BCL-1

10% SDS PAGE

GTXT111212 diluted at 1:3000



For full product information, images and publications, please visit our [website](http://www.genetex.com).