

Myelin basic protein antibody [1]

Cat. No. GTX11160

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Application	IHC-Fr, ELISA
Reactivity	Human, Rat, Bovine

Package
1 ml

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-Fr	1/10
ELISA	1/200-1/1,000

Not tested in other applications.

PROPERTIES

Form	Liquid
Buffer	Tissue culture supernatant, 0.2M Tris-HCl pH7.4, 5-10% foetal calf serum
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Bovine Myelin Basic Protein (MBP).
Purification	Tissue Culture Supernatant
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).