

RPS3A antibody

Cat. No. GTX111935

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Rat, Zebrafish

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 30 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.73 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human RPS3A Recombinant protein.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX111935 IHC-P Image

Immunohistochemical analysis of paraffin-embedded zebrafish tissue, using RPS3A antibody (GTX111935) at 1:300 dilution.



GTX111935 WB Image

Sample (30 µg of whole cell lysate)

A: zebrafish eye

12% SDS PAGE

GTX111935 diluted at 1:1000


GTX111935 WB Image


Sample (30 ug of whole cell lysate)

A: A431

B: Raji

12% SDS PAGE

GTX111935 diluted at 1:1000


GTX111935 IHC-P Image

Immunohistochemical analysis of paraffin-embedded U87 xenograft, using RPS3A antibody (GTX111935) at 1:500 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



For full product information, images and publications, please visit our [website](https://www.genetex.com).