

ZHX2 antibody [C1C3]

Cat. No. GTX112232

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-Fr, IP, ChIP-seq
Reactivity	Human, Mouse

References (10)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-Fr	Assay dependent
IP	1:100-1:500
ChIP-seq	Assay dependent

Not tested in other applications.

Calculated MW 92 kDa. ([Note](#))**Product Note** KO/KD validation is based on published data (PMID: 36037364).

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.7 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human ZHX2. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx112232 WB Image


Various whole cell extracts (30 µg) were separated by 5% SDS-PAGE, and the membrane was blotted with ZHX2 antibody [C1C3] (GTx112232) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.


GTx112232 IP Image

Immunoprecipitation of ZHX2 protein from MCF-7 whole cell extracts using 5 µg of ZHX2 antibody [C1C3] (GTx112232).

Western blot analysis was performed using ZHX2 antibody [C1C3] (GTx112232).

EasyBlot anti-Rabbit IgG (GTx221666-01) was used as a secondary reagent.


GTx112232 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using ZHX2(GTx112232) antibody at 1:200 dilution.



For full product information, images and publications, please visit our [website](https://www.genetex.com).