

Kinectin 1 antibody [C2C3], C-term

Cat. No. GTX112532

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

★★★★☆ Review (1)

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 156 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.34 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human Kinectin 1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx112532 IHC-P Image

Immunohistochemical analysis of paraffin-embedded Cal27 xenograft, using Kinectin 1A(GTx112532) antibody at 1:100 dilution.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx112532 WB Image

Sample (30 ug of whole cell lysate)

A: H1299

5% SDS PAGE

GTx112532 diluted at 1:1000



GTx112532 IHC-P Image

Kinectin 1A antibody [C2C3], C-term detects Kinectin 1A protein at cytoplasm in mouse brain by immunohistochemical analysis.

Sample: Paraffin-embedded mouse brain.

Kinectin 1A antibody [C2C3], C-term (GTx112532) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx112532 ICC/IF Image

Immunofluorescence analysis of methanol-fixed HeLa, using Kinectin 1A(GTx112532) antibody at 1:200 dilution.



For full product information, images and publications, please visit our [website](https://www.genetex.com).