

# Mad2L2 (3-14) peptide

**Cat. No. GTX11254**

**Species** Human

**Package**

50 µg

## Applications

**Product Note** This product can be used as the blocking peptide for GTX23630.

## Properties

**Form** Liquid

**Buffer** PBS

**Preservative** 0.02% Sodium azide

**Storage** Store vial at -20 °C prior to opening. Aliquot contents and freeze at -20 °C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use.

**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)

**Region/Sequence** Synthetic peptide, corresponding to amino acids 3-14 of Human Mad2L2.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).