

p23 antibody [N1C3-2]

Cat. No. GTX112655

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

References (4)

★★★★☆ Review (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 19 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.31 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human p23 Recombinant protein.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2

DATA IMAGES

GTX112655 WB Image

Sample (30 ug of whole cell lysate)


A: A431 (GTX27909)

B: H1299

C: Hela


12% SDS PAGE

GTX112655 diluted at 1:1000


GTX112655 IHC-P Image

Immunohistochemical analysis of paraffin-embedded HT-29 xenograft, using p23(GTX112655) antibody at 1:100 dilution.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX112655 WB Image

Sample (50 ug of whole cell lysate)

A: mouse brain

12% SDS PAGE

GTX112655 diluted at 1:1000


GTX112655 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed A431, using P23 (GTX112655) antibody at 1:200 dilution.



For full product information, images and publications, please visit our [website](https://www.genetex.com).