

beta Tubulin antibody

Cat. No. GTX112659

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 50 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human beta Tubulin. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx112659 ICC/IF Image

Immunofluorescence analysis of methanol-fixed HeLa, using beta Tubulin(GTx112659) antibody at 1:500 dilution.




GTx112659 WB Image

Sample (30 ug of whole cell lysate)
A: A549
10% SDS PAGE
GTx112659 diluted at 1:1000



GTx112659 IHC-P Image

Immunohistochemical analysis of paraffin-embedded HeLa xenograft, using beta Tubulin(GTx112659) antibody at 1:500 dilution.
Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx112659 WB Image

Sample (whole cell lysate)
A: 293T 20ug
B: 293T 10ug
C: 293T 5ug
10% SDS PAGE
GTx112659 diluted at 1:10000



For full product information, images and publications, please visit our [website](https://www.genetex.com).