

RAC2 antibody [N1C3]

Cat. No. GTX113262

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	Assay dependent
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 21 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.38 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human RAC2. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2


DATA IMAGES

GTx113262 IHC-P Image

RAC2 antibody [N1C3] detects RAC2 protein at cytoplasm by immunohistochemical analysis.

Sample: Paraffin-embedded mouse spleen.

RAC2 stained by RAC2 antibody [N1C3] (GTx113262) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTx113262 IHC-P Image

RAC2 antibody [N1C3] detects RAC2 protein at cytosol on mouse fore brain by immunohistochemical analysis.

Sample: Paraffin-embedded mouse fore brain.


RAC2 antibody [N1C3] (GTx113262) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTx113262 IHC-P Image

Immunohistochemical analysis of paraffin-embedded NCI-N87 xenograft, using RAC2(GTx113262) antibody at 1:500 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx113262 ICC/IF Image

RAC2 antibody [N1C3] detects RAC2 protein at cell membrane and cytoplasm by immunofluorescent analysis.

Sample: MOLT-4 cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: RAC2 stained by RAC2 antibody [N1C3] (GTx113262) diluted at 1:500.

Blue: Fluoroshield with DAPI (GTx30920).


GTx113262 WB Image

Various whole cell extracts (30 µg) were separated by 12% SDS-PAGE, and the membrane was blotted with RAC2 antibody [N1C3] (GTx113262) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).