

Sodium/Potassium ATPase beta 1 antibody

Cat. No. GTX113390

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

References (14)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 35 kDa. ([Note](#))

Product Note This antibody was raised against human Sodium Potassium ATPase Beta 1 Extracellular domain.

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the Extracellular domain of human Sodium Potassium ATPase Beta 1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 09 Page 1 of 2


DATA IMAGES

GTx113390 IHC-P Image

Sodium Potassium ATPase Beta 1 antibody detects ATP1B1 protein at membrane on human hepatoma by immunohistochemical analysis.

Sample: Paraffin-embedded hepatoma tissue.


Sodium Potassium ATPase Beta 1 antibody (GTx113390) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTx113390 IHC-P Image

Immunohistochemical analysis of paraffin-embedded MDAMB468 xenograft, using ATPase beta1(Na+/K+) (GTx113390) antibody at 1:100 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx113390 IHC-P Image

Sodium Potassium ATPase Beta 1 antibody detects ATP1B1 protein at membrane on Cal27 xenograft by immunohistochemical analysis.


Sample: Paraffin-embedded Cal27 xenograft.

Sodium Potassium ATPase Beta 1 antibody (GTx113390) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTx113390 WB Image

Mouse tissue extract (50 µg) was separated by 10% SDS-PAGE, and the membrane was blotted with Sodium/Potassium ATPase beta 1 antibody (GTx113390) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.


GTx113390 WB Image

Rat tissue extract (50 µg) was separated by 10% SDS-PAGE, and the membrane was blotted with Sodium/Potassium ATPase beta 1 antibody (GTx113390) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).