

hnRNP C1/C2 antibody

Cat. No. GTX113463

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Human, Mouse

References (6)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000
IP	Assay dependent

Not tested in other applications.

Calculated MW 34 kDa. ([Note](#))

Product Note KO/KD validation is based on published data (PMID: 31147722).

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 1% BSA, 10% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human hnRNP C1/C2. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES

GTX113463 IHC-P Image

hnRNP C1/C2 antibody detects hnRNP C1/C2 protein at nucleus in mouse brain by immunohistochemical analysis.

Sample: Paraffin-embedded mouse brain.

hnRNP C1/C2 antibody (GTX113463) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX113463 WB Image


Sample (30 µg of whole cell lysate)

A: Raji

12% SDS PAGE

GTX113463 diluted at 1:1000

The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.


GTX113463 WB Image

hnRNP C1/C2 antibody detects HNRNPC protein by western blot analysis.

A. 30 µg Neuro2A whole cell lysate/extract

B. 30 µg NIH-3T3 whole cell lysate/extract

C. 30 µg BCL-1 whole cell lysate/extract

D. 30 µg Raw264.7 whole cell lysate/extract

E. 30 µg C2C12 whole cell lysate/extract

12% SDS-PAGE

hnRNP C1/C2 antibody (GTX113463) dilution: 1:1000

The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).