

PSMA6 antibody

Cat. No. GTX113530

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 27 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human PSMA6. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2


DATA IMAGES

GTX113530 IHC-P Image

Proteasome 20S alpha 6 antibody detects Proteasome 20S alpha 6 protein at cytosol and nucleus on mouse lung by immunohistochemical analysis.


Sample: Paraffin-embedded mouse lung.

Proteasome 20S alpha 6 antibody (GTX113530) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTX113530 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed A431, using 20S Proteasome alpha6 (GTX113530) antibody at 1:200 dilution.


GTX113530 IHC-P Image

Immunohistochemical analysis of paraffin-embedded H1299 xenograft , using 20S Proteasome alpha6 (GTX113530) antibody at 1:100 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTX113530 WB Image

Proteasome 20S alpha 6 antibody detects PSMA6 protein by Western blot analysis.

A. 50 µg rat liver lysate/extract

12 % SDS-PAGE

Proteasome 20S alpha 6 antibody (GTX113530) dilution: 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).