

Factor IX antibody [N3C3]

Cat. No. GTX113689

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 52 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human Factor IX. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx113689 IHC-P Image

Factor IX antibody [N3C3] detects Factor IX protein at cytoplasm in human colon cancer by immunohistochemical analysis.

Sample: Paraffin-embedded human colon cancer.

Factor IX antibody [N3C3] (GTx113689) diluted at 1:250.



GTx113689 IHC-P Image

Factor IX antibody [N3C3] detects Factor IX protein at cytoplasm in mouse pancreas by immunohistochemical analysis.

Sample: Paraffin-embedded mouse pancreas.

Factor IX antibody [N3C3] (GTx113689) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTx113689 IHC-P Image

Factor IX antibody [N3C3] detects Factor IX protein at cytoplasm in rat liver by immunohistochemical analysis.

Sample: Paraffin-embedded rat liver.

Factor IX antibody [N3C3] (GTx113689) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx113689 IHC-P Image

Factor IX antibody [N3C3] detects Factor IX protein at cytoplasm by immunohistochemical analysis.

Sample: Paraffin-embedded mouse liver.

Factor IX stained by Factor IX antibody [N3C3] (GTx113689) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).