

TRIM32 antibody

Cat. No. GTX113937

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

References (5)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:20000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 72 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human TRIM32. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTx113937 IHC-P Image

TRIM32 antibody detects TRIM32 protein at cytoplasm in mouse brain by immunohistochemical analysis.

Sample: Paraffin-embedded mouse brain.

TRIM32 antibody (GTx113937) diluted at 1:400.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx113937 ICC/IF Image


TRIM32 antibody detects TRIM32 protein at cytoplasm by immunofluorescent analysis.

Sample: A375 cells were fixed in ice-cold MeOH for 5 min.

Green: TRIM32 protein stained by TRIM32 antibody (GTx113937) diluted at 1:500.

Blue: Hoechst 33342 staining.

Scale bar = 10 μ m.


GTx113937 WB Image


TRIM32 antibody detects TRIM32 protein by western blot analysis.

A. 30 μ g 293T whole cell lysate/extract

7.5% SDS-PAGE

TRIM32 antibody (GTx113937) dilution: 1:1000

The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.


GTx113937 WB Image

TRIM32 antibody detects TRIM32 protein by western blot analysis.

A. 5 μ g 293T whole cell lysate/extract

B. 5 μ g whole cell lysate/extract of DDDDK-human TRIM32-transfected 293T cells

7.5% SDS-PAGE

TRIM32 antibody (GTx113937) dilution: 1:20000

The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).