

Flightless I antibody [C2C3], C-term

Cat. No. GTX114410

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000

Not tested in other applications.

Calculated MW 145 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.95 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human Flightless I. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx114410 WB Image


Sample (30 ug of whole cell lysate)

A: Raji

5% SDS PAGE

Flightless I antibody

GTx114410 diluted at 1:1000




GTx114410 WB Image

Flightless I antibody [C2C3], C-term detects FLII protein by Western blot analysis.

A: 50 µg rat muscle lysate/extract

5 % SDS-PAGE

Flightless I antibody [C2C3], C-term (GTx114410) dilution: 1:1000




GTx114410 ICC/IF Image

Flightless I antibody [C2C3], C-term detects FLII protein at cytoplasm by immunofluorescent analysis.

Sample: Raji cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: FLII protein stained by Flightless I antibody [C2C3], C-term (GTx114410) diluted at 1:500.

Blue: Hoechst 33342 staining.



GTx114410 WB Image

Sample (50 ug of whole cell lysate)

A: Mouse brain

5% SDS PAGE

GTx114410 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).