

FADS1 antibody

Cat. No. GTX114528

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 52 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.91 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human FADS1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx114528 WB Image

Sample (30 ug of whole cell lysate)

A: A549

7.5% SDS PAGE

GTx114528 diluted at 1:1000



GTx114528 ICC/IF Image


FADS1 antibody detects FADS1 protein at mitochondria by immunofluorescent analysis.

Sample: Hep G2 cells were fixed in 2% paraformaldehyde/culture medium at 37°C for 30 min.

Green: FADS1 protein stained by FADS1 antibody (GTx114528) diluted at 1:500.

Blue: Hoechst 33342 staining.

Scale bar = 10 µm.



GTx114528 IHC-P Image

Immunohistochemical analysis of paraffin-embedded human colon carcinoma, using FADS1 (GTx114528) antibody at 1:250 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



For full product information, images and publications, please visit our [website](https://www.genetex.com).