

SMARCC2 / BAF170 antibody

Cat. No. GTX114778

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

 Review (1)Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.


Calculated MW 133 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.37 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the N-terminus region of human SMARCC2 / BAF170. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX114778 WB Image

Various whole cell extracts (30 µg) were separated by 5% SDS-PAGE, and the membrane was blotted with SMARCC2 antibody (GTX114778) diluted at 1:1000.




GTX114778 IHC-P Image

SMARCC2 antibody detects SMARCC2 protein at nucleolus in mouse duodenum by immunohistochemical analysis.

Sample: Paraffin-embedded mouse duodenum.

SMARCC2 antibody (GTX114778) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX114778 IHC-P Image

SMARCC2 antibody detects SMARCC2 protein at nucleolus in rat brain by immunohistochemical analysis.

Sample: Paraffin-embedded rat brain.

SMARCC2 antibody (GTX114778) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).