

ZNF133 antibody [N1C1-2]

Cat. No. GTX114896

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 73 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.54 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human ZNF133. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTx114896 IHC-P Image

ZNF133 antibody [N1C1-2] detects ZNF133 protein at nucleus in human breast cancer by immunohistochemical analysis.

Sample: Paraffin-embedded human breast cancer.

ZNF133 antibody [N1C1-2] (GTx114896) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTx114896 IHC-P Image

ZNF133 antibody [N1C1-2] detects ZNF133 protein at nucleus in human colon cancer by immunohistochemical analysis.


Sample: Paraffin-embedded human colon cancer.

ZNF133 antibody [N1C1-2] (GTx114896) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx114896 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using ZNF133(GTx114896) antibody at 1:200 dilution.


GTx114896 WB Image

Sample (30 ug of whole cell lysate)

A: HCT116 (GTx27909)

7.5% SDS PAGE

GTx114896 diluted at 1:5000



For full product information, images and publications, please visit our [website](http://www.genetex.com).