

ACOX3 antibody

Cat. No. GTX115077

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 78 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human ACOX3. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx115077 WB Image

Sample (30 ug of whole cell lysate)

A: A549

7.5% SDS PAGE

GTx115077 diluted at 1:1000




GTx115077 IHC-P Image

ACOX3 antibody detects ACOX3 protein at cytoplasm in human cervical carcinoma by immunohistochemical analysis.

Sample: Paraffin-embedded human cervical carcinoma.

ACOX3 antibody (GTx115077) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTx115077 WB Image

ACOX3 antibody detects ACOX3 protein by Western blot analysis.

A. 50 µg mouse liver lysate/extract

7.5 % SDS-PAGE

ACOX3 antibody (GTx115077) dilution: 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).