

SAFB2 antibody [N2C1], Internal

Cat. No. GTX115272

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 107 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.89 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human SAFB2. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX115272 ICC/IF Image

Immunofluorescence analysis of methanol-fixed A431, using Scaffold attachment factor B2(GTX115272) antibody at 1:200 dilution.


GTX115272 IHC-P Image

Immunohistochemical analysis of paraffin-embedded HBL435 xenograft, using Scaffold attachment factor B2(GTX115272) antibody at 1:100 dilution.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX115272 IHC-P Image

SAFB2 antibody [N2C1], Internal detects SAFB2 protein at nucleus in mouse duodenum by immunohistochemical analysis.

Sample: Paraffin-embedded mouse duodenum.

SAFB2 antibody [N2C1], Internal (GTX115272) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX115272 IHC-P Image

SAFB2 antibody [N2C1], Internal detects SAFB2 protein at nucleus in mouse brain by immunohistochemical analysis.

Sample: Paraffin-embedded mouse brain.

SAFB2 antibody [N2C1], Internal (GTX115272) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).