

U2AF65 antibody [C1C3]

Cat. No. GTX115622

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:5000-1:20000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 54 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.24 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human U2AF65. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX115622 IHC-P Image

U2AF65 antibody [C1C3] detects U2AF65 protein at cytoplasm and nucleus by immunohistochemical analysis.

Sample: Paraffin-embedded rat colon.

U2AF65 stained by U2AF65 antibody [C1C3] (GTX115622) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX115622 IHC-P Image

U2AF65 antibody [C1C3] detects U2AF65 protein at cytoplasm and nucleus by immunohistochemical analysis.

Sample: Paraffin-embedded mouse tissue.

U2AF65 stained by U2AF65 antibody [C1C3] (GTX115622) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX115622 ICC/IF Image

U2AF65 antibody [C1C3] detects U2AF65 protein at nucleus by immunofluorescent analysis.

Sample: A431 cells were fixed in 4% paraformaldehyde at RT for 15 min.


Green: U2AF65 stained by U2AF65 antibody [C1C3] (GTX115622) diluted at 1:500.

Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [GT114] (GTX628802) diluted at 1:1000.


GTX115622 IHC-P Image

Immunohistochemical analysis of paraffin-embedded human ovarian cancer, using U2AF65(GTX115622) antibody at 1:250 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTX115622 WB Image

Various whole cell extracts (30 µg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with U2AF65 antibody [C1C3] (GTX115622) diluted at 1:10000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).