

Thrombospondin antibody [TSP-B7]

Cat. No. GTX11579

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 129 kDa. ([Note](#))

Product Note

Specific for human platelet-derived thrombospondin and serum thrombospondin. Cross-reactivity has been observed with bovine thrombospondin in fetal calf serum (other species not tested). No cross-reactivity with human thrombin, fibrinogen, fibronectin, laminin or serum albumin is observed.

Properties

Form	Liquid
Buffer	Ascites
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Immunogen	Thrombin induced platelet releasate (Human).
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).