

Human Histone H2B protein, His tag


Cat. No. GTX115955-pro**Species** Human**Package**

100 µg, 25 µg

Properties

Form	Liquid
Buffer	PBS, 0.1% SDS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Region/Sequence	N-terminal His tagged full length protein of AAH96122
Expression System	E. coli
Purity	>90%
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES

**GTX115955-pro Image**

10µg of GTX115955-pro HIST1H2BD recombinant protein analyzed using SDS-PAGE and stained with coomassie blue and captured by black and white camera.



For full product information, images and publications, please visit our [website](#).