

PRMT5 antibody [N1C1]

Cat. No. GTX116004

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IP, ChIP assay
Reactivity	Human, Mouse

References (5)

★★★★★ Review (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000
IP	1:100-1:500
ChIP assay	Assay dependent

Not tested in other applications.

Calculated MW 73 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.31 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human PRMT5. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES

GTx116004 IHC-P Image

Immunohistochemical analysis of paraffin-embedded MCF7 xenograft, using PRMT5(GTx116004) antibody at 1:100 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTx116004 WB Image

Various whole cell extracts (30 µg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with PRMT5 antibody [N1C1] (GTx116004) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.



GTx116004 IHC-P Image

PRMT5 antibody [N1C1] detects PRMT5 protein at cytoplasm and nucleus by immunohistochemical analysis.

Sample: Paraffin-embedded mouse testis.

PRMT5 stained by PRMT5 antibody [N1C1] (GTx116004) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx116004 IP Image

PRMT5 antibody [N1C1] immunoprecipitates PRMT5 protein in IP experiments.

IP samples: Raji whole cell extract

A. 35 µg Raji whole cell extract

B. Control with 4 µg of preimmune Rabbit IgG

C. Immunoprecipitation of PRMT5 protein by 4 µg PRMT5 antibody [N1C1] (GTx116004)

7.5 % SDS-PAGE

The immunoprecipitated PRMT5 protein was detected by PRMT5 antibody [N1C1] (GTx116004) diluted at 1:500.

[EasyBlot anti-rabbit IgG (GTx221666-01) was used as a secondary reagent]



For full product information, images and publications, please visit our [website](https://www.genetex.com).