

Edc3 antibody [N2C2], Internal

Cat. No. GTX116386

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 56 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human Edc3. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx116386 WB Image


Edc3 antibody [N2C2], Internal detects EDC3 protein by Western blot analysis.

A. 30 µg NIH-3T3 whole cell lysate/extract

B. 30 µg C2C12 whole cell lysate/extract

7.5 % SDS-PAGE


Edc3 antibody [N2C2], Internal (GTx116386) dilution: 1:1000



GTx116386 IHC-P Image

Immunohistochemical analysis of paraffin-embedded MCF7 xenograft, using Edc3(GTx116386) antibody at 1:250 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min




GTx116386 IHC-P Image

Edc3 antibody [N2C2], Internal detects Edc3 protein at cytosol on mouse testis by immunohistochemical analysis.

Sample: Paraffin-embedded mouse testis.

Edc3 antibody [N2C2], Internal (GTx116386) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx116386 WB Image

Sample (30 ug of whole cell lysate)

A: A431 (GTx27909)

7.5% SDS PAGE

GTx116386 diluted at 1:3000



For full product information, images and publications, please visit our [website](https://www.genetex.com).