

SLC25A23 antibody

Cat. No. GTX116402

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 52 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human SLC25A23. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx116402 WB Image

Sample (50 ug of whole cell lysate)

A: mouse brain

10% SDS PAGE

GTx116402 diluted at 1:1000




GTx116402 ICC/IF Image

SLC25A23 antibody detects SLC25A23 protein at cytoplasm by immunofluorescent analysis.

Sample: SK-N-SH cells were fixed in 2% paraformaldehyde/culture medium at 37°C for 30 min.

Green: SLC25A23 protein stained by SLC25A23 antibody (GTx116402) diluted at 1:500.

Blue: Hoechst 33342 staining.




GTx116402 WB Image

Sample (30 ug of whole cell lysate)

A: IMR32

10% SDS PAGE

GTx116402 diluted at 1:500



GTx116402 IHC-P Image

Immunohistochemical analysis of paraffin-embedded human colon carcinoma, using SLC25A23(GTx116402) antibody at 1:250 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



For full product information, images and publications, please visit our [website](https://www.genetex.com).