

Human Hemoglobin Beta protein, His tag

Cat. No. GTX116412-pro

Species	Human
----------------	-------

Package

100 µg, 25 µg

Properties

Form	Liquid
-------------	--------

Buffer	PBS, 0.1% SDS
---------------	---------------

Preservative	No preservatives
---------------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
----------------	---

Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
----------------------	--

Region/Sequence	N-terminal His tagged full length protein of AAH07075
------------------------	---

Expression System	E. coli
--------------------------	---------

Purity	>90%
---------------	------

Conjugation	Unconjugated
--------------------	--------------

Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---



For full product information, images and publications, please visit our [website](#).