

MMP24 antibody

Cat. No. GTX11649

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-Fr, FCM, IP, ELISA
Reactivity	Human, Mouse

Package
25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2,000
IHC-Fr	Assay dependent
FCM	Assay dependent
IP	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 70 kDa. ([Note](#))

Product Note Reacts with reduced and native MMP-24. Recognizes only the pro-form of MMP-24 and not the active forms.

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	~1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	synthetic peptide corresponding to the propeptide region of mouse matrix metalloproteinase-24 (membrane-type matrix metalloproteinase-5).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).