

INTS4 antibody [N1N2], N-term

Cat. No. GTX116519

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 108 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.6 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the N-terminus region of human INTS4. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX116519 WB Image


Sample (30 ug of whole cell lysate)

A: U87-MG

B: MCF-7

7.5% SDS PAGE

GTX116519 diluted at 1:1000



GTX116519 IHC-P Image

INTS4 antibody [N1N2], N-term detects INTS4 protein at nucleus in mouse brain by immunohistochemical analysis.

Sample: Paraffin-embedded mouse brain.

INTS4 antibody [N1N2], N-term (GTX116519) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX116519 IHC-P Image

INTS4 antibody [N1N2], N-term detects INTS4 protein at nucleus in mouse ovary by immunohistochemical analysis.

Sample: Paraffin-embedded mouse ovary.

INTS4 antibody [N1N2], N-term (GTX116519) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX116519 IHC-P Image

Immunohistochemical analysis of paraffin-embedded AGS xenograft, using INTS4(GTX116519) antibody at 1:500 dilution.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](http://www.genetex.com).