

# PCDHA2 antibody [N1N2], N-term

**Cat. No. GTX116525**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Rabbit     |
| <b>Clonality</b>    | Polyclonal |
| <b>Isotype</b>      | IgG        |
| <b>Applications</b> | WB         |
| <b>Reactivity</b>   | Human      |

**Package**  
100 µl, 25 µl

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:1000-1:10000       |

Not tested in other applications.

**Calculated MW** 102 kDa. ([Note](#))


**Product Note** This antibody was raised against human PCDHA2 Extracellular domain.

## Properties

|                      |                                                                                                                                                                                                                                                                                                                                      |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                                                                                                                               |
| <b>Buffer</b>        | PBS, 1% BSA, 20% Glycerol                                                                                                                                                                                                                                                                                                            |
| <b>Preservative</b>  | 0.01% Thimerosal                                                                                                                                                                                                                                                                                                                     |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                           |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                             |
| <b>Immunogen</b>     | Recombinant protein encompassing a sequence within the Extracellular domain of human PCDHA2. The exact sequence is proprietary.                                                                                                                                                                                                      |
| <b>Purification</b>  | Purified by antigen-affinity chromatography.                                                                                                                                                                                                                                                                                         |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                                                                                                                                         |
| <b>Note</b>          | For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.<br><br>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product. |



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTx116525 WB Image**

Sample (30 ug of whole cell lysate)

A: NT2D1

B: IMR32

C: U87-MG

D: MCF-7

7.5% SDS PAGE

GTx116525 diluted at 1:5000



For full product information, images and publications, please visit our [website](https://www.genetex.com).