

Rat Anti-Mouse IgM (mu chain) antibody [LO-MM-3]

Cat. No. GTX11658

Host	Rat
Clonality	Monoclonal
Isotype	IgM
Application	ELISA
Reactivity	Mouse

Package
250 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	5µg/ml (as coating antibody)

Not tested in other applications.

Product Note

This antibody recognizes murine IgG3 (gamma heavy chain) and does not cross react with other murine immunoglobulin classes or subclasses.

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Polyclonal mouse IgM.
Purification	Purified by affinity chromatography From tissue culture supernatant
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).