

Visfatin antibody [N1N3]

Cat. No. GTX117444

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Rat

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 56 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human Visfatin. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx117444 WB Image

Visfatin antibody detects Visfatin protein by western blot analysis.

A. 30 μ g Jurkat whole cell lysate/extract


B. 30 μ g Raji whole cell lysate/extract

C. 30 μ g K562 whole cell lysate/extract

D. 30 μ g NCI-H929 whole cell lysate/extract

7.5 % SDS-PAGE

Visfatin antibody (GTx117444) dilution: 1:1000



GTx117444 WB Image


Sample (30 ug of whole cell lysate)

A: A431

B: Raji

7.5% SDS PAGE

GTx117444 diluted at 1:1000



GTx117444 WB Image

Visfatin antibody detects Visfatin protein by western blot analysis.


A. 30 μ g A549 whole cell lysate/extract

B. 30 μ g HCT116 whole cell lysate/extract

C. 30 μ g MCF-7 whole cell lysate/extract

7.5 % SDS-PAGE

Visfatin antibody (GTx117444) dilution: 1:1000



GTx117444 IHC-P Image

Visfatin antibody [N1N3] detects Visfatin protein at cytoplasm in rat liver by immunohistochemical analysis.

Sample: Paraffin-embedded rat liver.

Visfatin antibody [N1N3] (GTx117444) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](http://www.genetex.com).