

CRIP2 antibody [N1C3-2]

Cat. No. GTX117712

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 22 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.55 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human CRIP2. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX117712 ICC/IF Image

Immunofluorescence analysis of methanol-fixed HeLa, using CRIP2(GTX117712) antibody at 1:200 dilution.




GTX117712 IHC-P Image

CRIP2 antibody [N1C3-2] detects CRIP2 protein at cytoplasm on human lung adenocarcinoma blood vessel by immunohistochemical analysis.

Sample: Paraffin-embedded lung adenocarcinoma blood vessel.

CRIP2 antibody [N1C3-2] (GTX117712) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTX117712 WB Image

Sample (30 ug of whole cell lysate)

A: HeLa

12% SDS PAGE

GTX117712 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).