

# COPS4 antibody [N3C3]

**Cat. No. GTX117852**

|                     |                   |
|---------------------|-------------------|
| <b>Host</b>         | Rabbit            |
| <b>Clonality</b>    | Polyclonal        |
| <b>Isotype</b>      | IgG               |
| <b>Applications</b> | WB, ICC/IF, IHC-P |
| <b>Reactivity</b>   | Human             |

**Package**  
100 µl, 25 µl

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500-1:3000         |
| ICC/IF             | 1:100-1:1000         |
| IHC-P              | 1:100-1:1000         |

Not tested in other applications.

**Calculated MW** 46 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS, 1% BSA, 20% Glycerol                                                                                                                                                                                                  |
| <b>Preservative</b>  | 0.01% Thimerosal                                                                                                                                                                                                           |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| <b>Immunogen</b>     | Recombinant protein encompassing a sequence within the center region of human COPS4. The exact sequence is proprietary.                                                                                                    |
| <b>Purification</b>  | Purified by antigen-affinity chromatography.                                                                                                                                                                               |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTx117852 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using COPS4(GTx117852) antibody at 1:200 dilution.



### GTx117852 WB Image

Sample (30 ug of whole cell lysate)  
A: HCT116  
10% SDS PAGE  
GTx117852 diluted at 1:1000



### GTx117852 IHC-P Image

Immunohistochemical analysis of paraffin-embedded MDA-MB157 xenograft, using COPS4(GTx117852) antibody at 1:100 dilution.  
Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).