

MMP27 antibody

Cat. No. GTX11788

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
 25 µg

Applications
Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1,000

Not tested in other applications.

Calculated MW 59 kDa. ([Note](#))

Properties
Form Liquid

Buffer PBS, 50% Glycerol

Preservative 0.05% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration ~1? mg/ml (Please refer to the vial label for the specific concentration.)

Immunogen synthetic peptide corresponding to the hinge region of human matrix metalloproteinase-27.

Purification Purified by affinity chromatography

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


 For full product information, images and publications, please visit our [website](#).