

ASB8 antibody

Cat. No. GTX117973

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 32 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human ASB8. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX117973 WB Image

Sample (30 ug of whole cell lysate)

A: IMR32

12% SDS PAGE

GTX117973 diluted at 1:1000




GTX117973 IHC-P Image

ASB8 antibody detects ASB8 protein at cytoplasm in rat skin by immunohistochemical analysis.

Sample: Paraffin-embedded rat skin.


ASB8 antibody (GTX117973) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX117973 ICC/IF Image

Immunofluorescence analysis of methanol-fixed HeLa, using ASB8(GTX117973) antibody at 1:500 dilution.



GTX117973 IHC-P Image

Immunohistochemical analysis of paraffin-embedded U87 xenograft, using ASB8(GTX117973) antibody at 1:500 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



For full product information, images and publications, please visit our [website](https://www.genetex.com).