

RNF212 antibody

Cat. No. GTX118251

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 33 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human RNF212. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX118251 IHC-P Image

RNF212 antibody detects RNF212 protein at cytoplasm in human esophagus cancer by immunohistochemical analysis.

Sample: Paraffin-embedded human esophagus cancer.

RNF212 antibody (GTX118251) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX118251 IHC-P Image

RNF212 antibody detects RNF212 protein at cytoplasm in human breast cancer by immunohistochemical analysis.

Sample: Paraffin-embedded human breast cancer.


RNF212 antibody (GTX118251) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX118251 IHC-P Image

Immunohistochemical analysis of paraffin-embedded U87 xenograft, using RNF212(GTX118251) antibody at 1:500 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTX118251 WB Image

Sample (30 ug of whole cell lysate)

A: PC-3

12% SDS PAGE

GTX118251 diluted at 1:1000



For full product information, images and publications, please visit our [website](http://www.genetex.com).